

# Originate™

## The Origin of Great Restorations.

Color coded for easy identification of size, upper, and lower

Palate-free to reduce the "gag reflex" for more patient comfort and much easier model work.


Patented material locking plate holds in impression material.

Rigid design insures no side-to-side movement associated with bite trays resulting in more accurate impressions.

Rough inner surface for material retention

Perforations allow the impression material to flow into the tray and lock in.

Sturdy, ergonomic handle.



**IT-101**  
Upper Small



**IT-102**  
Upper Medium



**IT-103**  
Upper Large



**IT-201**  
Lower Small



**IT-202**  
Lower Medium



**IT-203**  
Lower Large

# Originate™

## Great Restorations Begin With Great Impressions!

You have upgraded your impression materials over the years. Now it's time to upgrade your impression tray.

Originate™ Disposable Impression Trays are palate-free which aids in reducing the gag reflex. These trays are full arch with enough rigidity to prevent the side-to-side movement associated with flimsy bite trays. Originate™ can be used with any type of impression material used in fabricating crown and bridge restorations, bleaching trays and mouthguards.

Impression material flows into the unique locking plate on the bottom of the tray locking it securely into the tray.

Each Originate™ tray is individually packaged and manufactured with Microban® Antimicrobial Product Protection to increase sterility. Trays are color-coded to easily indicate size and to identify upper and lower trays. Each tray also includes an identification tag.

You have upgraded your impression materials over the years. Now it is time to upgrade your impression tray. Make the switch to Originate™ Disposable Impression Trays today.



**Originate™ Assortment Pack**  
Order No. IT-ASST

**The Originate™ Assortment Pack is a great way to start!**

Each assortment pack includes (6) IT-101 Upper Small, (6) IT-102 Upper Medium, (2) IT-103 Upper Large, (2) IT-201 Lower Small, (6) IT-202 Lower Medium, and (2) IT-203 Lower Large.

Contact an AXIS Dental Product Specialist or your preferred authorized dealer representative today.



US Patent 5,890,895 and 6,045,359. Patent Pending (Canada).

**Product Specialist:**  
**888.654.2947**  
www.axisdental.com  
L0507-AXORIG